
AGIERRE[®]

Engineering & Innovative Solutions

*Advanced Technologies
in the Handling of Powder*

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Vacuum Transfer System for powders and grains

AGIERRE[®]
Engineering & Innovative Solutions

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The Agierre's Pneumatic Transfer Systems are the result of the most advanced technology in the handling of powders.

They are mainly applied both in the productive department (on mixers, dryers, fluid bed coolers, blending systems) and in the packaging units (on packaging machines, tableting press machines, etc).

The functioning is based on the suction (by vacuum) of powders and grains that must be transferred from a container to another one.

The Vacuum Pump (modelled according to the Venturi Principle) creates vacuum and depression inside the Vacuum Conveyor and consequently the product suction. When the product is inside the Vacuum Conveyor, it will be separated from the air by a filter. After that, the product falls down thanks to gravity and centrifugal force. The main functioning is based on the compressed air, so the contamination of products is minimum. The Vacuum allows the handling of the powder with a low velocity which guaranties a constant and consistent handling of powders without shock and de-mixing.

On the Vacuum Transfer System can be installed other pumps such us liquid ring pump or similar.

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The Vacuum Transfer Systems have an innovative design and they can be customized according to the technical requirements of each Customer. According to the technical and processing specification the Vacuum Transfer Systems can have different configurations:

- Basic;
- Overpressure;
- Inertisation;

The connections could be with clamp, wing nuts or flanged (also in the UNI PN16 version).

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*AUTOMATIC DISPENSING PLANT
FOR RAW MATERIALS*

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The Vacuum Conveyors (patent pending) are manufactured in compliance with the "cGMP guidelines", the UNI EN ISO norms and the Atex and Machine Directives. According to the AGIERRE's manufacturing philosophy, all Conveyors can be customized in according to specific needs.

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The Conveyor could be also used in the applications with explosive powders which require, during the discharge phase, an oxygen residual value inferior to 6%. It's possible with the inertisation process.

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Vacuum Transfer System with C.I.P. and S.I.P. systems

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These cleaning systems find their main application in the pharmaceutical and chemical industries, particularly in the sterile production areas where it's needed to preserve the chemical and physical properties of the product and avoid the "cross contamination".

The Vacuum Transfer Systems can be equipped with some accessories that allow to wash, dry and sterilize it through the C.I.P. (cleaning in place) and S.I.P. (steam in place) systems.

The C.I.P. system cleans the filter and all internal surfaces of Vacuum Conveyor by rotating nozzles (patent pending system).

The C.I.P. system is based on a mechanical cleaning system that ensure the complete removal of all product residues inside the Vacuum Conveyor. After washing phase all surfaces can be dried with air and/or nitrogen. This Conveyor has an innovative and specific design without any point of product or liquid accumulation.

The Vacuum Conveyors can be managed by a pneumatic control console or by a PLC system.

They are manufactured with stainless steel AISI 316L and are in compliance with the "cGMP guidelines", the UNI EN ISO norms, the Atex and Machine Directives.

According to the Agierre's manufacturing philosophy the VacuumConveyors with C.I.P. system can be customized according to the specific needs of each Customer.

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Vacuum Transfer System in FILTERVAC version

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Tank for powder

The Vacuum Conveyor type FILTERVAC is used to charge machines like Mixer, Reactor and Bicone where it's not possible to install the Vacuum Conveyor directly on the machine. The main functioning is based on the product suction by vacuum created inside the machine.

The Pump installed on the FILTERVAC generates the vacuum in the Machine that allows the continuous product suction from a drum, bin or other machines.

During the charging phase the FILTERVAC makes also the filter cleaning. At the end of the suction process is made the so called "purge phase" consisting in the powder recovery that can be done manually or automatically through the suction hopper installed under the discharge point of the FILTERVAC.

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Through the Ferrule of Vapour Collector the Vacuum Conveyor can be installed directly on reactors and/or dissolvers even with both solvents and vapour at 100 °, thanks to a structure with double cavity wall which creates a barrier of laminar flow. The double cavity wall traps the vapour of reactor which goes up to the inlet connection of reactor. This system avoids any product accumulation on the discharge valve of the Vacuum Conveyor installed on the reactor, simplifying the discharging phase of product.

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The Lump Breaker is manufactured in accordance with the Machinery Directive 2006/42/EC and upon request it can be manufactured for installation in ATEX zones.

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The Lump Breaker model AGR-BREAKER allows to crush the lumps that could form in the materials during the production process or the handling of powders. The Lump Breaker Systems are particularly suitable for hygroscopic materials or those that tend to clog, and thanks to their versatility can be used in all industrial sectors.

It is characterized by a self-cleaning design with very small dimensions and can also be mounted within existing systems or systems for emptying (bags, drums and Big-Bag), or usually located immediately above vessel, filter dryers, centrifuges, mixers and generally on the process machine.

The breaking of lumps takes place through the rotation of special blades within a grid of contrast moved by means of a gearmotor.

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The lump breaker aspiration lance is designed to work in vacuum and to be connected with our Vacuum Conveyor. It finds the main application in the chemical and pharmaceutical Industries for the handling of humid powders with big lumps that the standard suction lance can't transport.

The lump breaker lance can be manufactured also with a fluidisation air/nitrogen system which guarantee the right air/nitrogen proportion for the fluid dynamics balance.

According to the Agierre's manufacturing philosophy the lump breaker lance is manufactured in compliance with the "cGMP guidelines", UNI EN ISO and safety norms.

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Vacuum Transfer System PED (Pressure Rated)

AGIERRE®

Engineering & Innovative Solutions

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Vacuum Transfer System Pressure Rated

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All Vacuum Conveyors PED are manufactured according to the "cGMP guidelines", the UNI EN ISO norms, ATEX Directive, Machinery Directive and PED Directive 97/23/EC.

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The PED Conveyors are manufactured for the need to install the Vacuum Transfer System on reactors PED classified (according to the Directive 97/23/EC).

The main specifications are:

- Material : AISI 316L (Stainless Steel);
- Thickness of plate: 30/10;
- Internal and external finish: mirror polished;
- Connection between body and lid: flange, wing nut or clamp;
- Functioning: pneumatic , electropneumatic or with PLC
- Pressure: -1+3 bargs / -1+6 bargs

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
Also the PED Conveyors can be manufactured with C.I.P. (Cleaning in place) and S.I.P. (Steam in place) Systems. This type of Vacuum Conveyor is named GMA.

- | | |
|---|--|
| <p>1 Connessione Pompa e ingresso prodotto : Valvola a farfalla
<i>Connection Vacuum Pump and inlet connection : Butterfly valve</i></p> <p>2 Connessione aria/azoto in pressione: Valvola a Membrana/Farfalla
<i>Connection air/nitrogen: Diaphragm/Butterfly Valve</i></p> <p>3 Valvola di Sicurezza
<i>Safety Valve</i></p> | <p>4 Filtro Acciaio Inox
<i>Stainless Steel Filter</i></p> <p>5 Filtro in Acciaio inox con sistema di lavaggio per C.I.P. e S.I.P.
<i>Stainless Steel Filter with the system of cleaning for C.I.P. and S.I.P.</i></p> <p>6 Valvola a Farfalla sanitaria per ambienti sterili
<i>Sanitary Butterfly Valve for Sterile Areas</i></p> |
|---|--|

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According to the AGIERRE's projecting and manufacturing philosophy, also this type of Vacuum Conveyor (PED) can be customized according to the needs of the Customer.

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Vacuum Transfer System in HOPPERVAC version

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The Vacuum Conveyor type HOPPERVAC is mainly used when it's necessary to transport the product into a tank and/or container using the HOPPERVAC as a storage unit to supply another machine. In this case the HOPPERVAC is used not only as a powder transfer system but also as a container of product accumulation.

The functioning is the same of standard Vacuum Conveyor based on the product suction by the vacuum pump.

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The HOPPERVAC like all other AGIERRE's Vacuum Transfer Systems can be manufactured with the C.I.P. system (Cleaning in Place) which allows the complete cleaning of the filter and the internal walls of Vacuum Conveyor.

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According to the AGIERRE's projecting philophy, the HOPPERVAC can be customized according to the needs of customer, always in compliance with the "cGMP guidelines", the UNI EN ISO norm, the Machine and ATEX Directives and with the safety norm.

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
***HOPPERVAC ro
charge mill and
dosing system***

Applicazioni Conveyor in Vacuum *Vacuum Conveyor Installations*

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The Lifting Column for Pharma BIN is mainly used in the chemical, pharmaceutical and food industries to optimize the production process of the handling of powder and avoid the "cross-contamination" processes. Upon request it can be manufactured also according to the Atex Directive therefore it can be installed in powdery zones. It has a primary internal galvanized iron structure and an external stainless steel AISI 304 carter at which is connected the arm.

The Lifting system consists of a "worm screw" moved by a self-breaking gear motor. All functions of the Lifting Column (horizontal and vertical movement as well as all customizable technical specifications) are managed automatically through the control panel. The lifting system is adjusted by superior, intermediate and inferior limit switches that can be adjusted according to the specific application and process requirements. The rotation (horizontal movement) may be manual or automatically managed by the control console.

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The Lifting Column for BIN has a specific arm for the installation of BIN. The Control Panel with PLC system manages automatically the Vibrator activation and the opening and closing system of bottom valve.

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The Lifting Columns for BIN and the BIN Blender Columns can be manufactured in ATEX version too.

Carrello di sostegno Bin
Truck to support Bin

The BIN Blender Column can be also used as a Lifting column; it has a stainless steel AISI 304 external structure and a galvanized carbon steel internal and primary structure.

The Lifting system is made of a "worm screw" and "wheel gear" while the rotation system is made by a self-braking epicyclic gear motor which guaranty power and elegancy.

The lifting system is adjusted by limit switches which can adjusted the maximum and minimum lifting height of the arm according to the specific application and process requirements. The arm has a specific connection for Bin with 10° angle for a better mixing process. The forks of Bin Blender Column have an automatic device of lock/unlock of BIN (patent pending). All functions of the BIN Blender Column (horizontal and vertical movement as well as all customizable technical specifications) are managed automatically by the control panel with PLC system.

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Lifting Column for Vacuum Transfer Systems

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The Lifting Columns of AGR-LTFT series are mainly used to lift and support Vacuum Conveyor. They have the same structure and functioning of the Standard Lifting Columns of AGR-LIFT series. All functions of the Lifting Column (horizontal and vertical movement as well as all customizable technical specifications) are managed automatically by the control panel with PLC system.

The lifting system is adjusted by the superior, intermediate and inferior limit switches that allow the right and perfect position of the Vacuum Conveyor on the Machine and the charge point. They are mainly used in the following application and installation such as:

- where there is the need to feed Machines where the point of load is high;
- when the Vacuum Conveyor can be not installed on the Machine;
- To load Schet Filling Machines, Horizontal Ribbon Mixer, Dryers, etc;
- For safety reasons to prevent the operator dismantle the Vacuum Conveyor installed at elevated heights.

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The Lifting Columns for Vacuum Conveyor can be also mobile.

They are mainly used when:

- it is not possible to do a fixed installation;
- the same machine can be used on several and different applications in different production areas.

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The Lifting Columns for Drums are mainly used to lift and rotate carton and plastic drums. They have the same structure and functioning of the Standard Lifting Columns of AGR-LIFT series.

Also in this type of Lifting Column the lifting system is made of a "worm screw".

All functions of the Lifting Column for Drums (horizontal and vertical movement) are managed automatically by the control panel that could have a pneumatic, electrical-pneumatic or PLC system.

The Lifting Columns for Drums are used for the charge and feed of machine where is needed to do the handling of powders or tablets by gravity.

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All Lifting Columns can be manufactured also in ATEX version

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The discharge of product is made without any powders emission through a closing and tightening system properly designed for the specific Big-Bag.

The whole system is controlled by PLC system and is designed in accordance with the Machinery Directive 2006/42/EC and all the safety systems; upon request it is designed for applications in ATEX zones.

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Big Bag emptying structure with
load system with Lift Truck

The Big-Bag emptying structure is a modular system with different configurations, depending on the type of application, used for the empty Big-Bag system.

The structure can be made of painted steel or satin AISI 304 and consists of two main parts: the base of fixed structure to the ground and a structure that supports the Big-Bag.

On the top of the structure can have two types of configuration for the Big-Bag support: the structure for Big-Bag using the forklift and an alternative structure to support the beam, where there is an electric hoist motor with a crane to support the Big-Bag. The Big-Bag Emptying structure has pneumatic movable supports on which place the Big-Bag; during the suction phase of the product, the pneumatic actuators activate the movable supports with a rotation movement with alternate torque that allows the rapid, accurate, and complete emptying of the Big-Bag.

The system can be placed on load cells by which is possible to see the actual weight contained in the Big-Bag and the weight of the product transferred and/or dosed in the receiving system.

Colonna Svuota Big Bag

Lifting Column for emptying Big Bag

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The Lifting Column for Big Bag has a movable emptying system with bag-wrench that ensures the full powder containment thanks to the closing system with two closing handles.

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The Lifting Column for emptying Big Bag of AGR-LIFT series is designed and manufactured for the handling and supporting Big Bag.

It is mainly used in the Chemical, Pharmaceutical and Food Industry especially in those processes of handling of powder to optimize the production flow and avoid cases of "cross pollution". Upon request it could be installed in powdery environments and in ATEX classification zone.

The column has a stainless steel AISI 304 external structure and a primary galvanized carbon steel internal structure at which is fixed the arm which supports the load.

All functions of the lifting column (horizontal and vertical movement as well as all its customizable specification) are automatically managed by the control panel installed on the Column. The rotation and lifting systems is managed by a superior, intermediate and inferior limit switches adjustable according to the specific installation and process requirements.

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*The Tumbler is designed according to the cGMP standards norms and in accordance with the current safety regulations.
It is provided with a set of magnetic sensors and photoelectric barriers that guarantee both the passive and active safety simultaneously.*

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*The Tumbler is designed with an high-tech design and a structure that makes it easy to clean and to look inside it.
It consists of:*

- **Supporting structure:** it consists of carbon steel profiles anchored to the floor by means of anchor bolts.
- **Lifting platform:** it consists of a stainless steel AISI 304 arm connected to the sliding guides moved by the screw-nut screw system driven by an electric motor.
- **External Carter:** it consists of a stainless steel AISI 304 carter shaped to cover the supporting structure of the Tumbler.
- **Gear motor:** it consists of an electric motor coupled to a reduction gear for the rotation function of the bin.
- **Control panel**

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The Bin Blender is designed with a highly technological design and a structure that makes it easy to make maintenance activities and clean. It is made of the following components:

- *Galvanized carbon steel internal structure that represents the machine carrier structure. The arm which support load and the gearmotor are fixed on this structure.*
- *External stainless steel AISI 304 Casing.*
- *Stainless steel AISI 304 security Frame with lexan panels which encompass the bin to be mixed. The entry in the mixing area is allowed by two doors equipped with electro-locking system.*
- *Anti-static wheels with stainless steel AISI 304 support that allow the easy movement of the machine within the department; The rear casters are fitted with a brake.*
- *Stainless steel AISI 304 Arm designed to blend Bin of 5 liters to 300 liters, with a maximum capacity of 150 kgs.*
- *Gearmotor group.*
- *Touch panel from which it's possible to carry out the loading/unloading and mixing of the bin.*

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The Bin Blender is designed in full compliance with current "cGMP" standard norms and in compliance with current safety regulations, UNI EN ISO standards and ATEX regulations (when required).

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The manual bag emptying system represents the ideal solution for emptying of powder or granules in the spaces with the complete absence of dust. It consists of a grid with a support front shelf mounted on a hopper supported by four feet and a hopper with door protection. The system can be realized with different types of material with a high degree of finish in the version with or without filter. It can be managed by a control panel with pneumatic operating system or with the PLC according to the specific application.

The system is so designed to minimize the internal stagnation of the product and can be customized with different solutions to meet the wide multiplicity and variety of applications.

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Opening system of door: the system can be customized with an openable glass door with the application of two gloves that allows to manipulate the product internally without the environment or the operator are exposed to the product.

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Depending on the type of application on the discharge point of the hopper can be also installed a valve or a breaker both connected to a small hopper connected to a powder transport system.

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


Ventilation system: the upper chamber can be connected to a general extraction system or can be used without ventilation only by means of the installed filters.

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
Filtration system: the working chamber is provided with a special filtration system according to the degree of containment that is required (polyester filters or HEPA filters).

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Lifting and tilting drums system: according to the types of application the system can be arranged for emptying of drums thanks to an automatic tilting mechanism for drums with an inflatable seal ensuring ab OEB4 level of containment;

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Continuous liner system: it can be installed on a side of the upper chamber and used for the safe removal of waste or for the introduction of instruments into the chamber. Inside there is an hinged door or a cap to guarantee the closure of such connection when the system is not used.

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Emptying Systems for bags and drums

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Manual emptying system for bags connected to a Vacuum transfer system

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Automatic emptying system for drums

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Finitura interna lucidata a specchio
Internal finish: mirror polished

Carrello amovibile
Movable truck

BIN in cGMP execution

Size

- from 100 to 2.000 litres.

Finish

- Internal : mirror polished $Ra < 0,1\mu$ with the removal of weldings;
- External : satin-finished $Ra < 0,6\mu$ with the removal of weldings;

Materials

- Stainless Steel (AISI 316L) for the internal surfaces.
- Stainless Steel (AISI 304) for the external surfaces.

Standard Connection of the discharge point

- Clamp connection

Accessories

- TC connection on the discharge point.
- Pneumatic Vibrator to make easier the discharge of powder.

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
Accessory:

- *Pneumatic Vibrator with Suction Cup*

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Bottom Valve with clamp connection.

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
External finish: satin-finished with the weldings removal and mirror polished

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On request our BIN can be installed on wheels

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
Closing Tri-Clamp for Lid with safety Pin.

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Lid with a stainless steel Filter and TC hatch for the connection to the other machineries.

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
Accessories for Lid:

- Ferrule
- Ferrule for the discharge of product

Sistema svuotamento fusti con trasporto polveri in ambiente classificato OEB4/5
Drums Emptying System with Powder Transfer System in areas OEB4/5 classified

Applicazioni Applications

Sistema integrato di Colonna Alza Bin e Gira Bin per carico mulino
Integrated system of Bin Blender Column for the charge of mill

Struttura svota Big Bag con trasporto polveri in ambiente confinato OEB4/5
Emptying Big Bag Structure with Powder Transfer System in OEB4/5 areas