



SafeLock™

Automatic needle guard for
pre-filled syringe

Volume 

Viscosity 

Device Type

Primary Container

Dosage Type

Fill Volume

Route of Administration

Handling

Injection Feedback

Disposable Safety Syringe

1mL Long Syringe with Staked Needle or
2.25mL Syringe with Staked Needle

Fixed | Full Contents

0.1-1mL | 1-2.25mL

Subcutaneous

Manual, one-handed use, automatic lock-out at end of dose

Viewing window

Handling Instructions



Remove cap



Press & hold against injection site vertically



Press down plunger rod to administer the drug, causing needle cover to lock into place



Discard in sharps container

Key Features

- ✓ Simple 3-step operation
- ✓ Visual feedback
- ✓ One-handed operation
- ✓ Compact design
- ✓ Locking needle cover
- ✓ Large labeling area
- ✓ Preconfigured, with colors customized to client specifications
- ✓ Compatible with Rigid Needle Shield